



OPA1 Antibody, Biotin conjugated

Product Code	CSB-PA016340HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60313
Immunogen	Recombinant Human Dynamin-like 120 kDa protein, mitochondrial protein (183-418AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dynamin-like 120 kDa protein, mitochondrial (EC 3.6.5.5) (Optic atrophy protein 1) [Cleaved into: Dynamin-like 120 kDa protein, form S1], OPA1, KIAA0567
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OPA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.