



PREP Antibody, Biotin conjugated

Product Code	CSB-PA018664LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48147
Immunogen	Recombinant Human Prolyl endopeptidase protein (539-691AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prolyl endopeptidase (PE) (EC 3.4.21.26) (Post-proline cleaving enzyme), PREP, PEP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PREP