



UNC13B Antibody, FITC conjugated

Product Code	CSB-PA025622LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14795
Immunogen	Recombinant Human Protein unc-13 homolog B protein (12-215AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein unc-13 homolog B (Munc13-2) (munc13), UNC13B, UNC13
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	UNC13B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.