



POU4F1 Antibody, Biotin conjugated

Product Code	CSB-PA22749D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01851
Immunogen	Recombinant Human POU domain, class 4, transcription factor 1 protein (258-419AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	POU domain, class 4, transcription factor 1 (Brain-specific homeobox/POU domain protein 3A) (Brain-3A) (Brn-3A) (Homeobox/POU domain protein RDC-1) (Oct-T1), POU4F1, BRN3A RDC1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	POU4F1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.