



ITIH1 Antibody, Biotin conjugated

Product Code	CSB-PA4779D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19827
Immunogen	Recombinant Human Inter-alpha-trypsin inhibitor heavy chain H1 protein (403-638AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inter-alpha-trypsin inhibitor heavy chain H1 (ITI heavy chain H1) (ITI-HC1) (Inter-alpha-inhibitor heavy chain 1) (Inter-alpha-trypsin inhibitor complex component III) (Serum-derived hyaluronan-associated protein) (SHAP), ITIH1, IGHEP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ITIH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.