



MTMR3 Antibody, Biotin conjugated

Product Code	CSB-PA622675HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13615
Immunogen	Recombinant Human Myotubularin-related protein 3 protein (652-899AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myotubularin-related protein 3 (EC 3.1.3.48) (FYVE domain-containing dual specificity protein phosphatase 1) (FYVE-DSP1) (Phosphatidylinositol-3,5-bisphosphate 3-phosphatase) (EC 3.1.3.95) (Phosphatidylinositol-3-phosphate phosphatase) (EC 3.1.3.64) (Zinc finger FYVE domain-containing protein 10), MTMR3, KIAA0371 ZFYVE10
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MTMR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.