



KPNB1 Antibody, HRP conjugated

Product Code	CSB-PA622929LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14974
Immunogen	Recombinant Human Importin subunit beta-1 protein (9-214AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Importin subunit beta-1 (Importin-90) (Karyopherin subunit beta-1) (Nuclear factor p97) (Pore targeting complex 97 kDa subunit) (PTAC97), KPNB1, NTF97
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	KPNB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.