



PKN3 Antibody, Biotin conjugated

Product Code	CSB-PA744059LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P5Z2
Immunogen	Recombinant Human Serine/threonine-protein kinase N3 protein (1-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase N3 (EC 2.7.11.13) (Protein kinase PKN-beta) (Protein-kinase C-related kinase 3), PKN3, PKNBETA
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PKN3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.