



# TNFRSF25 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA839000LB01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q93038
<b>Immunogen</b>	Recombinant Human Tumor necrosis factor receptor superfamily member 25 protein (221-417AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Tumor necrosis factor receptor superfamily member 25 (Apo-3) (Apoptosis-inducing receptor AIR) (Apoptosis-mediating receptor DR3) (Apoptosis-mediating receptor TRAMP) (Death receptor 3) (Lymphocyte-associated receptor of death) (LARD) (Protein WSL) (Protein WSL-1), TNFRSF25, APO3 DDR3 DR3 TNFRSF12 WSL WSL1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	TNFRSF25