



TNFRSF25 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA839000LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q93038 |
| Immunogen | Recombinant Human Tumor necrosis factor receptor superfamily member 25 protein (221-417AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Tumor necrosis factor receptor superfamily member 25 (Apo-3) (Apoptosis-inducing receptor AIR) (Apoptosis-mediating receptor DR3) (Apoptosis-mediating receptor TRAMP) (Death receptor 3) (Lymphocyte-associated receptor of death) (LARD) (Protein WSL) (Protein WSL-1), TNFRSF25, APO3 DDR3 DR3 TNFRSF12 WSL WSL1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | TNFRSF25 |