




CELSR3 Antibody

Product Code	CSB-PA852817LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYQ7
Immunogen	Recombinant Human Cadherin EGF LAG seven-pass G-type receptor 3 protein (514-718AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:200-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cadherin EGF LAG seven-pass G-type receptor 3 (Cadherin family member 11) (Epidermal growth factor-like protein 1) (EGF-like protein 1) (Flamingo homolog 1) (hFmi1) (Multiple epidermal growth factor-like domains protein 2) (Multiple EGF-like domains protein 2), CELSR3, CDHF11 EGFL1 FMI1 KIAA0812 MEGF2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CELSR3

Image



IHC image of CSB-PA852817LA01HU diluted at 1:300 and staining in paraffin-embedded human liver cancer performed on a Leica Bond™ system. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4°C overnight. The primary is detected by a biotinylated secondary antibody and visualized using an HRP conjugated SP system.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.