



CELSR3 Antibody, FITC conjugated

Product Code	CSB-PA852817LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYQ7
Immunogen	Recombinant Human Cadherin EGF LAG seven-pass G-type receptor 3 protein (514-718AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cadherin EGF LAG seven-pass G-type receptor 3 (Cadherin family member 11) (Epidermal growth factor-like protein 1) (EGF-like protein 1) (Flamingo homolog 1) (hFmi1) (Multiple epidermal growth factor-like domains protein 2) (Multiple EGF-like domains protein 2), CELSR3, CDHF11 EGFL1 FMI1 KIAA0812 MEGF2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CELSR3