



CELSR3 Antibody, Biotin conjugated

Product Code	CSB-PA852817LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYQ7
Immunogen	Recombinant Human Cadherin EGF LAG seven-pass G-type receptor 3 protein (514-718AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cadherin EGF LAG seven-pass G-type receptor 3 (Cadherin family member 11) (Epidermal growth factor-like protein 1) (EGF-like protein 1) (Flamingo homolog 1) (hFmi1) (Multiple epidermal growth factor-like domains protein 2) (Multiple EGF-like domains protein 2), CELSR3, CDHF11 EGFL1 FMI1 KIAA0812 MEGF2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CELSR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.