



MYO18A Antibody, HRP conjugated

Product Code	CSB-PA856414LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92614
Immunogen	Recombinant Human Unconventional myosin-XVIIIa protein (1849-2054AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Unconventional myosin-XVIIIa (Molecule associated with JAK3 N-terminus) (MAJN) (Myosin containing a PDZ domain) (Surfactant protein receptor SP-R210) (SP-R210), MYO18A, CD245 KIAA0216 MYSPDZ
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYO18A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.