



RUBCN Antibody, Biotin conjugated

Product Code	CSB-PA856415LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92622
Immunogen	Recombinant Human Run domain Beclin-1-interacting and cysteine-rich domain-containing protein (518-733AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Run domain Beclin-1-interacting and cysteine-rich domain-containing protein (Rubicon) (Beclin-1 associated RUN domain containing protein) (Baron), RUBCN, KIAA0226
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RUBCN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.