



NSMAF Antibody, HRP conjugated

Product Code	CSB-PA856416LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92636
Immunogen	Recombinant Human Protein FAN protein (1-233AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAN (Factor associated with neutral sphingomyelinase activation) (Factor associated with N-SMase activation), NSMAF, FAN
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NSMAF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.