



SYNGAP1 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA857006LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96PV0 |
| Immunogen | Recombinant Human Ras/Rap GTPase-activating protein SynGAP protein (1161-1343AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Ras/Rap GTPase-activating protein SynGAP (Neuronal RasGAP) (Synaptic Ras GTPase-activating protein 1) (Synaptic Ras-GAP 1), SYNGAP1, KIAA1938 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | SYNGAP1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |