



PAFAH2 Antibody, HRP conjugated

Product Code	CSB-PA857448LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99487
Immunogen	Recombinant Human Platelet-activating factor acetylhydrolase 2, cytoplasmic protein (1-206AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Platelet-activating factor acetylhydrolase 2, cytoplasmic (EC 3.1.1.47) (Serine-dependent phospholipase A2) (SD-PLA2) (hSD-PLA2), PAFAH2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PAFAH2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.