



# P2RX7 Antibody, Biotin conjugated

Product Code	CSB-PA859103LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99572
Immunogen	Recombinant Human P2X purinoceptor 7 protein (110-247AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	P2X purinoceptor 7 (P2X7) (ATP receptor) (P2Z receptor) (Purinergic receptor), P2RX7
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	P2RX7