



P2RX7 Antibody, Biotin conjugated

Product Code	CSB-PA859103LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99572
Immunogen	Recombinant Human P2X purinoceptor 7 protein (110-247AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	P2X purinoceptor 7 (P2X7) (ATP receptor) (P2Z receptor) (Purinerbic receptor), P2RX7
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	P2RX7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.