



LARS Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA873736LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9P2J5 |
| Immunogen | Recombinant Human Leucine--tRNA ligase, cytoplasmic protein (614-892AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Leucine--tRNA ligase, cytoplasmic (EC 6.1.1.4) (Leucyl-tRNA synthetase) (LeuRS), LARS, KIAA1352 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | LARS |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |