



ANKRD2 Antibody, Biotin conjugated

Product Code	CSB-PA875635LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZV1
Immunogen	Recombinant Human Ankyrin repeat domain-containing protein 2 protein (1-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ankyrin repeat domain-containing protein 2 (Skeletal muscle ankyrin repeat protein) (hArpp), ANKRD2, ARPP
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ANKRD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.