



WDR35 Antibody, HRP conjugated

Product Code	CSB-PA879071LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2L0
Immunogen	Recombinant Human WD repeat-containing protein 35 protein (954-1181AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat-containing protein 35 (Intraflagellar transport protein 121 homolog), WDR35, IFT121 KIAA1336
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WDR35
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.