



WDR35 Antibody, Biotin conjugated

Product Code	CSB-PA879071LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2L0
Immunogen	Recombinant Human WD repeat-containing protein 35 protein (954-1181AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat-containing protein 35 (Intraflagellar transport protein 121 homolog), WDR35, IFT121 KIAA1336
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WDR35