



VPS18 Antibody, FITC conjugated

Product Code	CSB-PA882173LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P253
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 18 homolog protein (298-515AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 18 homolog (hVPS18), VPS18, KIAA1475
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VPS18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.