



TPT1P8 Antibody, HRP conjugated

Product Code	CSB-PA884505LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAU6
Immunogen	Recombinant Human Putative translationally-controlled tumor protein-like protein TPT1P8 protein (1-139AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative translationally-controlled tumor protein-like protein TPT1P8 (Putative apoptosis inhibitor FKSG2), TPT1P8
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TPT1P8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.