



TPT1P8 Antibody, Biotin conjugated

Product Code	CSB-PA884505LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAU6
Immunogen	Recombinant Human Putative translationally-controlled tumor protein-like protein TPT1P8 protein (1-139AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative translationally-controlled tumor protein-like protein TPT1P8 (Putative apoptosis inhibitor FKSG2), TPT1P8
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TPT1P8