



ADNP Antibody, Biotin conjugated

Product Code	CSB-PA887964HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2P0
Immunogen	Recombinant Human Activity-dependent neuroprotector homeobox protein (867-1102AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Activity-dependent neuroprotector homeobox protein (Activity-dependent neuroprotective protein), ADNP, ADNP1 KIAA0784
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ADNP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.