



RFNG Antibody, HRP conjugated

| | |
|---------------------|--|
| Product Code | CSB-PA896922LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9Y644 |
| Immunogen | Recombinant Human Beta-1,3-N-acetylglucosaminyltransferase radical fringe protein (71-331AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Beta-1,3-N-acetylglucosaminyltransferase radical fringe (EC 2.4.1.222) (O-fucosylpeptide 3-beta-N-acetylglucosaminyltransferase), RFNG |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | RFNG |