



ASAH2B Antibody

Product Code	CSB-PA002171ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C7U1
Immunogen	Recombinant Homo sapiens ASAH2B protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative inactive neutral ceramidase B (ASAH2-like protein) (Putative inactive N-acylsphingosine amidohydrolase 2B) (Putative inactive non-lysosomal ceramidase B), ASAH2B, ASAH2C ASAH2L
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ASAH2B
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.