



Camk4 Antibody

Product Code	CSB-PA004471ZA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13234
Immunogen	Recombinant Rattus norvegicus Camk4 protein
Raised In	Rabbit
Species Reactivity	Rattus norvegicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium/calmodulin-dependent protein kinase type IV (CaMK IV) (EC 2.7.11.17) (CaM kinase-GR) (Calspermin), Camk4
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	Camk4
Product Modify	Unmodified