



Capn1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA004490ZA01MO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O35350 |
| Immunogen | Recombinant Mus musculus Capn1 protein |
| Raised In | Rabbit |
| Species Reactivity | Mus musculus |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Calpain-1 catalytic subunit (EC 3.4.22.52) (Calcium-activated neutral proteinase 1) (CANP 1) (Calpain mu-type) (Calpain-1 large subunit) (Micromolar-calpain) (muCANP), Capn1, Canp1 Capa1 |
| Immunogen Species | Mus musculus (Mouse) |
| Lead Time | Made-to-order (12-14 weeks) |
| Gene Names | Capn1 |
| Product Modify | Unmodified |