



CD8A Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA004966ZA01BO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P31783 |
| Immunogen | Recombinant Bovine T-cell surface glycoprotein CD8 alpha chain (CD8A), partial (26-189AA) |
| Raised In | Rabbit |
| Species Reactivity | Bos taurus |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Protein A/G |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | T-cell surface glycoprotein CD8 alpha chain (CD antigen CD8a), CD8A |
| Immunogen Species | Bos taurus (Bovine) |
| Lead Time | Made-to-order |
| Gene Names | CD8A |
| Product Modify | Unmodified |