



CHRNA1 Antibody

Product Code	CSB-PA005386ZA01TOT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02710
Immunogen	Recombinant Tetronarce californica CHRNA1 protein
Raised In	Rabbit
Species Reactivity	Tetronarce californica
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Tetronarce californica (Pacific electric ray) (Torpedo californica)
Lead Time	Made-to-order (12-14 weeks)
Target Names	CHRNA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.