



# DUSP14 Antibody

<b>Product Code</b>	CSB-PA007243ZA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O95147
<b>Immunogen</b>	Recombinant Homo sapiens DUSP14 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Homo sapiens
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Dual specificity protein phosphatase 14 (EC 3.1.3.16) (EC 3.1.3.48) (MKP-1-like protein tyrosine phosphatase) (MKP-L) (Mitogen-activated protein kinase phosphatase 6) (MAP kinase phosphatase 6) (MKP-6), DUSP14, MKP6
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Gene Names</b>	DUSP14
<b>Product Modify</b>	Unmodified