



EPO Antibody

Product Code	CSB-PA007743ZA01HO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q867B1
Immunogen	Recombinant Equus caballus EPO protein
Raised In	Rabbit
Species Reactivity	Equus caballus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Erythropoietin, EPO
Immunogen Species	Equus caballus (Horse)
Lead Time	Made-to-order (12-14 weeks)
Target Names	EPO
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.