



# GHRH Antibody

|                     |   |
|---------------------|---|
| Product Code        | CSB-PA009412ZA01SH  |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| Uniprot No.         | P07217  |
| Immunogen           | Recombinant Ovis aries GHRH protein   |
| Raised In           | Rabbit  |
| Species Reactivity  | Ovis aries  |
| Tested Applications | ELISA, WB (ensure identification of antigen)  |
| Form                | Liquid  |
| Conjugate           | Non-conjugated  |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                      |
| Purification Method | Protein A/G   |
| Isotype             | IgG   |
| Clonality           | Polyclonal  |
| Alias               | Somatoliberin (Growth hormone-releasing factor) (GRF) (Growth hormone-releasing hormone) (GHRH), GHRH |
| Immunogen Species   | Ovis aries (Sheep)  |
| Lead Time           | Made-to-order (12-14 weeks)   |
| Gene Names          | GHRH  |
| Product Modify      | Unmodified  |