



GNPDA1 Antibody

Product Code	CSB-PA009630ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46926
Immunogen	Recombinant Homo sapiens GNPDA1 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucosamine-6-phosphate isomerase 1 (EC 3.5.99.6) (Glucosamine-6-phosphate deaminase 1) (GNPDA 1) (GlcN6P deaminase 1) (Oscillin), GNPDA1, GNPI HLN KIAA0060
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Target Names	GNPDA1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.