



Gsta1 Antibody

Product Code	CSB-PA009970ZA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00502
Immunogen	Recombinant Rattus norvegicus Gsta1 protein
Raised In	Rabbit
Species Reactivity	Rattus norvegicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutathione S-transferase alpha-1 (EC 2.5.1.18) (GST 1-1) (GST 1a-1a) (GST A1-1) (GST B) (Glutathione S-transferase Ya-1) (GST Ya1) (Ligandin), Gsta1
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Gsta1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.