



HSP90AA1 Antibody

Product Code	CSB-PA010802ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07900
Immunogen	Recombinant Homo sapiens HSP90AA1 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock protein HSP 90-alpha (Heat shock 86 kDa) (HSP 86) (HSP86) (Lipopolysaccharide-associated protein 2) (LAP-2) (LPS-associated protein 2) (Renal carcinoma antigen NY-REN-38), HSP90AA1, HSP90A HSPC1 HSPCA
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	HSP90AA1
Product Modify	Unmodified