



IFNG Antibody

Product Code	CSB-PA011050ZA01SH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17773
Immunogen	Recombinant Ovis aries IFNG protein
Raised In	Rabbit
Species Reactivity	Ovis aries
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon gamma (IFN-gamma), IFNG
Immunogen Species	Ovis aries (Sheep)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	IFNG
Product Modify	Unmodified