



IL6 Antibody

Product Code	CSB-PA011664ZA01SH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29455
Immunogen	Recombinant Ovis aries IL6 protein
Raised In	Rabbit
Species Reactivity	Ovis aries
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-6 (IL-6), IL6
Immunogen Species	Ovis aries (Sheep)
Lead Time	Made-to-order (12-14 weeks)
Target Names	IL6
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.