



# LHB Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA012910ZA01BO  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | P04651  |
| <b>Immunogen</b>           | Recombinant Bos taurus LHB protein  |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Bos taurus  |
| <b>Tested Applications</b> | ELISA, WB (ensure identification of antigen)  |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Non-conjugated  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                              |
| <b>Purification Method</b> | Antigen Affinity Purified   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Lutropin subunit beta (Luteinizing hormone subunit beta) (LH-B) (LSH-B) (LSH-beta) (Lutropin beta chain), LHB |
| <b>Immunogen Species</b>   | Bos taurus (Bovine)   |
| <b>Lead Time</b>           | Made-to-order (12-14 weeks)   |
| <b>Gene Names</b>          | LHB   |
| <b>Product Modify</b>      | Unmodified  |