



MAPT Antibody

Product Code	CSB-PA013481ZA01MOW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57786
Immunogen	Recombinant Macaca mulatta MAPT protein
Raised In	Rabbit
Species Reactivity	Macaca mulatta
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Microtubule-associated protein tau (Neurofibrillary tangle protein) (Paired helical filament-tau) (PHF-tau), MAPT, TAU
Immunogen Species	Macaca mulatta (Rhesus macaque)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	MAPT
Product Modify	Unmodified