



PABPC1 Antibody

Product Code	CSB-PA017352ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11940
Immunogen	Recombinant Homo sapiens PABPC1 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polyadenylate-binding protein 1 (PABP-1) (Poly(A)-binding protein 1), PABPC1, PAB1 PABP1 PABPC2
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	PABPC1
Product Modify	Unmodified