



SAA1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA020656ZA01HO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P19857 |
| Immunogen | Recombinant Equus caballus SAA1 protein |
| Raised In | Rabbit |
| Species Reactivity | Equus caballus |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Serum amyloid A protein (SAA) [Cleaved into: Amyloid protein A (Amyloid fibril protein AA)], SAA1 |
| Immunogen Species | Equus caballus (Horse) |
| Lead Time | Made-to-order (12-14 weeks) |
| Gene Names | SAA1 |
| Product Modify | Unmodified |