



TGFA Antibody

Product Code	CSB-PA023445ZA01SH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98135
Immunogen	Recombinant Ovis aries TGFA protein
Raised In	Rabbit
Species Reactivity	Ovis aries
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Protransforming growth factor alpha [Cleaved into: Transforming growth factor alpha (TGF-alpha) (EGF-like TGF) (ETGF) (TGF type 1)] (Fragment), TGFA, TGF-A
Immunogen Species	Ovis aries (Sheep)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	TGFA
Product Modify	Unmodified