



# TNFRSF1B Antibody

<b>Product Code</b>	CSB-PA023978ZA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P20333
<b>Immunogen</b>	Recombinant Homo sapiens TNFRSF1B protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Homo sapiens
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Tumor necrosis factor receptor superfamily member 1B (Tumor necrosis factor receptor 2) (TNF-R2) (Tumor necrosis factor receptor type II) (TNF-RII) (TNFR-II) (p75) (p80 TNF-alpha receptor) (CD antigen CD120b) (Etanercept) [Cleaved into: Tumor necrosis factor receptor superfamily member 1b, membrane form; Tumor necrosis factor-binding protein 2 (TBP-2) (TBPII)], TNFRSF1B, TNFBR TNFR2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Target Names</b>	TNFRSF1B
<b>Product Modify</b>	Unmodified
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.