



TXNDC17 Antibody

Product Code	CSB-PA025376ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRA2
Immunogen	Recombinant Homo sapiens TXNDC17 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Thioredoxin domain-containing protein 17 (14 kDa thioredoxin-related protein) (TRP14) (Protein 42-9-9) (Thioredoxin-like protein 5), TXNDC17, TXNL5
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	TXNDC17
Product Modify	Unmodified