



Ier Antibody

Product Code	CSB-PA188874ZA01Ba
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O87544
Immunogen	Recombinant Lysobacter enzymogenes Ier protein
Raised In	Rabbit
Species Reactivity	Lysobacter enzymogenes
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoproteinase Arg-C (Ier)
Immunogen Species	Lysobacter enzymogenes
Lead Time	Made-to-order (12-14 weeks)
Gene Names	Ier
Product Modify	Unmodified